

Water covers over 70% of the Earth's surface...
and constitutes over 70% of a human body...



- Water Phosphorylcholine or PC has a strong affinity for water.

NEW! from
GRACE MEDICAL

- PC coating develops a water layer on it's surface when exposed to fluids.

- The water molecules bind themselves tightly around the ventilation tube making it difficult for proteins and cells to displace the water and bind to the material of the tube.

CAT. #	DESCRIPTION	I.D. LUMEN
510-001C	Donaldson Type, w/tab	1.14mm
510-021C	Shepard Grommet Type, w/tab	1.02mm
510-101C	T-Tube, 12mm L	1.14mm
510-111C	Modified T-Tube, 6mm L	1.27mm
510-121C	Triune	1.35mm
510-131C	Collar Button	1.27mm
510-141C	Taper Vent	1.14mm
510-051C	Paparella Type	1.02mm



GRACE MEDICAL

8500 Wolf Lake Drive, Suite 110 Memphis, TN 38133 • USA
Phone: 901.386.0990 • Toll Free 1.866.472.2363 • Fax: 901.386.0950
info@gracemedical.com • www.gracemedical.com