

# Dornhoffer



## Dornhoffer

- Low profile contoured HA head
- Malleus groove for additional stability
- Windows for increased visibility
- Footplate shoe can be adjusted to fit crus remnants
- Footplate shoe helps center and stabilize TORPs

REF. #	DESCRIPTION
635	ALTO Total
695	ALTO Partial
690	2mm Fixed Length Partial
637	ALTO Total w/ Windows
687	ALTO Partial w/ Windows
691	2mm Fixed Length Partial w/ Windows
626	ALTO Total w/ Windows and Integrated Shoe
636-080*	.8mm Shoe

\*The Footplate Shoe can be used with all titanium shaft total implants

Designed in conjunction with John Dornhoffer, MD, Little Rock, AR



Grace Medical We Listen.

©2018 Grace Medical, Inc. All Rights Reserved.

8500 Wolf Lake Drive, Suite 110  
Memphis, TN 38133 USA  
Phone 1.901.386.0990  
Fax 1.901.386.0950  
Toll Free 1.866.472.2363  
info@gracemedical.com  
www.gracemedical.com

LIT0064 CID6029 Rev. 2018-09