



Silicone Cold Mold Technology



Grace Medical Cold Mold Silicone Tube Surface



Traditional Hot Mold Silicone Tube Surface



Cold Silicone Mold Technology

Cold Silicone → Cold Mold → Hot Oven

- Creates Smooth Tube Surface
- Provides Less Area For Bacteria Adhesion

Hot Silicone Mold Technology

Cold Silicone → Hot Mold → Hot Oven

- Creates “Orange Peel” Tube Surface
- Rougher Surface May Create Platform For Bacteria Adhesion



Cold Mold Technology is available on all Grace Medical silicone vent tubes.