

# Wildcat Four 2 Three Prosthesis

Windows for visibility



Centered shaft design

The Wildcat Four2Three enables conversion of a type 4 tympanoplasty to a type 3. The one piece design allows for easy insertion in a canal wall down with an absent stapes superstructure.

Hollow shaft



Solid shaft

Option of solid or hollow shaft



Cat.#	Description
707-300	Wildcat Fixed Length 3.00mm
707-350	Wildcat Fixed Length 3.50 mm
707-400	Wildcat Fixed Length 4.00 mm
707-450	Wildcat Fixed Length 4.50 mm
707-500	Wildcat Fixed Length 5.00 mm
707-550	Wildcat Fixed Length 5.50 mm
706-300	Wildcat Fixed Length 3.00mm (hollow shaft)
706-350	Wildcat Fixed Length 3.50 mm (hollow shaft)
706-400	Wildcat Fixed Length 4.00 mm (hollow shaft)
706-450	Wildcat Fixed Length 4.50 mm (hollow shaft)
706-500	Wildcat Fixed Length 5.00 mm (hollow shaft)
706-550	Wildcat Fixed Length 5.50 mm (hollow shaft)

Designed in conjunction with Alan Micco, MD, Chicago, IL, Cliff Megerian, MD, Cleveland, OH, Seilesh Babu, Farmington Hills, MI  
8500WolfLakeDrive,Suite110•Memphis,TN38133•USATollFree1.866.472.2363•www.gracemedical.com