

2in1 ALTO



Grace Medical We Listen.

2in1 ALTO

- Total and Partial in one prosthesis
- TORP functional length 3.00mm to 7.00mm
- PORP functional length 1.25mm to 5.75mm
- Modified base to accommodate both stapes footplate and crura
- Multidirectional head for optimal contact with TM

REF. #	DESCRIPTION
642	2in1 ALTO



Grace Medical We Listen.

©2020 Grace Medical, Inc. All Rights Reserved.

8500 Wolf Lake Drive, Suite 110
Memphis, TN 38133 USA
Phone 1.901.386.0990
Fax 1.901.386.0950
Toll Free 1.866.472.2363
info@gracemedical.com
www.gracemedical.com

LIT0241 CID7047 Rev. 2020-08