


Wildcat Pro ALTO Total



Grace Medical We Listen.

Wildcat Pro ALTO Total

- Adjustable length from 3mm to 7mm
- Multidirectional head for optimal contact with TM
- Large base for increased stability
- Crura cutouts to accommodate crura remnants

REF. #	DESCRIPTION
613	Wildcat Pro ALTO Total



Grace Medical We Listen.

©2020 Grace Medical, Inc. All Rights Reserved.

8500 Wolf Lake Drive, Suite 110
Memphis, TN 38133 USA
Phone 1.901.386.0990
Fax 1.901.386.0950
Toll Free 1.866.472.2363
info@gracemedical.com
www.gracemedical.com

LIT0242 CID7048 Rev. 2020-08